

VELOCITY[®]
orthopedics

**Value
Through
Innovation**

Clear 'N Dry[™]

Velocity[®]...

A premium line of
instrumentation for
arthroscopic surgery.



*Set includes disposable cannula
and disposable conical obturator.*

**Control Flow, Stay Dry,
See Clearly**

Control Flow, Stay Dry, See Clearly

Clear 'N Dry™

- Latex free seal - control flow and stay dry
- Transparent cannula for optimal monitoring of instrumentation and suture
- Unique threaded design for easy insertion and secure retention
- Luer and 1/4" tube connection on an adjustable fluid port
- Color coding by diameter
- Packaged sterile, single-use, ten per box
- Set includes disposable cannula and disposable conical obturator



ORDERING INFORMATION

3100	6.5mm x 72mm	Arthroscopy Cannula w/ stopcock
3101	8.5mm x 72mm	Arthroscopy Cannula w/ stopcock



Velocity® Orthopedics, Inc.

10700 Jersey Blvd., Suite 360, Rancho Cucamonga, CA 91730 USA Tel: 909.987.4343 Fax: 909.987.5800

© 2017 Velocity Orthopedics, Inc.

velocityorthopedics.com